




NARISHIGE SELECTION CHART FOR HUMAN IVF SYSTEM

Microscope	Micromanipulator Set	Adapter	Microinjector
IX3 Special Set for Evident	ON4MD, <u>Motorized Set for both sides,</u> E0437550, NAR-MMAN-MOT-IX3-ON4MD (without adapter) 2x 3-axis Hanging Joystick Oil Hydraulic Micromanipulators, 2x Motorized Manipulators 	ON-73MD, E0437551, NAR-AD-ON-73MD 	<u>YOU NEED ONE FOR EACH SIDE</u> IM-11-2MD, <u>Pneumatic Microinjector,</u> E0437558, NAR-MINJ-AIR-IM-11-2MD 
	MTK-1MD, <u>Set for both sides (Takanome),</u> E0437552, NAR-MMAN-IX-MTK-1MD (without adapter) 2x 4-axis Hanging Joystick Oil Hydraulic Micromanipulators, 2x Hydraulic manipulators 	IX2: NK-O2MD, E0437553, NAR-AD-MTK-IX2-NK-O2MD  IX3: NK-O3MD, E0437554, NAR-AD-MTK-IX3-NK-O3MD 	<u>OR/AND</u> IM-21MD, <u>Oil Microinjector,</u> E0437559 NAR-MINJ-OIL-IM-21MD 
	MOR-1MD, <u>Motorized Set for both sides,</u> E0437555, NAR-MMAN-MOT-IX-MOR-1MD (without adapter) 2x 3-axis Hanging Joystick Oil Hydraulic Micromanipulator, 2x Motorized Manipulators 	IX2: NO-PIX-3MD, E0437556, NAR-AD-MOR-IX2-NO-PIX-3MD  IX3: NO-PIX-5MD, E0437557, NAR-AD-MOR-IX3-NO-PIX-5MD 	NOTE! You can use instead of Narishige the microinjectors from Calibre (Eppendorf): Celltram 4m Air or Celtram 4m Oil (E0437181, EP-CELLTRAM-4M-AIR OR E0437182, EP-CELLTRAM-4M-OIL)

Additional products
SL-ICSI SYSTEM
(spindle observation)
for semi-motorized IX73 or IX83 configured for IVF



HEATING PLATE
Okolab or Tokai Hit



CAMERA



ANTIVIBRATION PROTECTION
Table Newport
(E0437072, NP-ACAFTA-7501200HG)



LASERS
(Available in some countries. Please contact PMM.)

Complete system from Narishige:

Option 1 (Special Set for Evident, Motorized):
Micromanipulator ON4MD,
Adapter for IX3: ON-73MD,
2x Microinjectors:
IM-11-2MD and/or IM-21MD



Option 2 (Takanome, Hydraulic):
Micromanipulator MTK-1MD,
Adapter for IX2: NK-O2MD
OR for IX3: NK-O3MD,
2x Microinjectors:
IM-11-2MD and/or IM-21MD



Option 3 (Motorized):
Micromanipulator MOR-1MD,
Adapter for IX2: NO-PIX-3MD
OR for IX3: NO-PIX-5MD,
2x Microinjectors:
IM-11-2MD and/or IM-21MD

